1521 E Third Street, Charlotte, NC 28204 | 1-800-633-3914 | howco.com



ITEM# 58001 HOWCO[™] DYNABRITE[™]2.0

| # OF DROPS | RATIO | SUGGESTED USAGE RANGE |
|------------|--------|-----------------------|
| 55 | 30 : 1 | Tough Cleaning** |
| 50 | 33 : 1 | Tough Cleaning** |
| 45 | 37 : 1 | Normal Cleaning |
| 40 | 41:1 | Normal Cleaning |
| 35 | 47 : 1 | Easy Cleaning |

At the above dilutions this product offers > 20 Grains per Gallon (gpg) of water hardness control.

**In seasonal winter conditions, use tough cleaning titration levels to maintain a high degree of cleaning.

1521 E Third Street, Charlotte, NC 28204 | 1-800-633-3914 | howco.com



Card

titration

ITEM# 58001 **HOWCO[™] DYNABRITE ™ 2.0**

 Rinse vial three times with sample to be tested.
Fill vial to 10 mL.



- 2. Add 3 drops Phenolphthalein (R-0638) and swirl to mix. The sample should turn PINK.
- Add 3 drops of Total Alkalinity (R-0645) and swirl to mix. The sample should turn PURPLE.

- 4. Add Sulfuric Acid (R-0736) 1 drop at a time. BEGIN COUNTING NUMBER OF DROPS. While holding the bottle vertically, swirl between each drop. The color will change from PURPLE to GREEN.
- 5. Continue adding the Sulfuric Acid (R-0736) 1 drop at a time. Swirl between each drop. Stop when the color changes from GREEN to PINK.



NOTE the number of drops and refer to the titration card for recommendations.